

TORQ[®] NP



Conexão
Digital Implant

INNOVATIVE DESIGN FOR FAST AND PRECISE INSERTION, WITH TECHNOLOGY
THAT ENSURES BONE COMPACTION AND FORCE DISSIPATION

Short - 977004

Long - 977005



CHARACTERISTICS

General Characteristics

- Reduced collar to increase bone thickness at the implant's cervical area
- Threads with high cutting power
- Internal angulation: 24°
- Internal thread: M 1.6
- Web[®] (enhanced bone contact)
- Bioplatform[®] (switch platform type)

Main Indications

- Single, partial, and multiple implantations
- Conventional protocol and immediate loading
- Dual indication: maxilla and mandible

VULCANO[®] Surface Treatment

Ø 3.5	Ø 3.75	Ø 4.0	Ø 5.0
545385 - 3,5x8,5	545785 - 3,75x8,5	545285 - 4,0x8,5	545585 - 5,0x8,5
545310 - 3,5x10	545710 - 3,75x10	545210 - 4,0x10	545510 - 5,0x10
545311 - 3,5x11,5	545711 - 3,75x11,5	545211 - 4,0x11,5	545511 - 5,0x11,5
545313 - 3,5x13	545713 - 3,75x13	545213 - 4,0x13	545513 - 5,0x13
545315 - 3,5x15	545715 - 3,75x15		



Recommended Insertion
Torque **60N.cm**



Hard[®]

Titânio Hard[®] 955(MPa)
(73% stronger)

POROUS[®] Surface Treatment

Ø 3.5	Ø 3.75	Ø 4.0	Ø 5.0
525385 - 3,5x8,5	525785 - 3,75x8,5	525285 - 4,0x8,5	525585 - 5,0x8,5
525310 - 3,5x10	525710 - 3,75x10	525210 - 4,0x10	525510 - 5,0x10
525311 - 3,5x11,5	525711 - 3,75x11,5	525211 - 4,0x11,5	525511 - 5,0x11,5
525313 - 3,5x13	525713 - 3,75x13	525213 - 4,0x13	525513 - 5,0x13
525315 - 3,5x15	525715 - 3,75x15		



INDEXED MORSE 24° - FLEX GOLD® NP - Ø3,5 - 3,75 - 4,3 - 5,0

| SHORT® NP - Ø3,75 - 4,3


FLASH® NP - Ø3,5 - 4,3 - 5,0

| TORQ® NP - Ø3,5 - 3,75 - 4,0 - 5,0

exocad®

exocad® / inLab

Probe
Gingival



Healing Caps



Abutment

Conventional
Transfer
(plaster)

Scan Conect
Digital
(Peek Medical)

DG
Analogues
(Digital / Plaster)

Fastener

Provisional

Cylinders
(copings)

Protective
Covers

Multiple			Unitary - Multiple		Unitary	Unitary - Multiple		Unitary				
MB - Multi Base			MU - Micro Unit Mini Pilar		EC - Esteticone	UMD - Universal Modular Digital		UB - UCLA Base CrCo		TiBase Small Large		
Φ3,8 Torque 30Ncm 230020 - 0,8mm 230021 - 1,5mm 230022 - 2,5mm 230023 - 3,5mm 230024 - 4,5mm			Φ4,8 Torque 30Ncm Solid 136211 - 1,5mm 136212 - 2,5mm 136213 - 3,5mm 132095 - 5,5mm 136123 - 0,5mm E 2 pieces 136101 - 1,5mm 136102 - 2,5mm 136103 - 3,5mm 17° 148102 - 2,5mm 148103 - 3,5mm 30° Linha 0 w/o Hex 032103 - 3,5mm 032104 - 4,5mm		Φ4,8 Torque 30Ncm 157104 022191 - 1,5mm 022192 - 2,5mm 022193 - 3,5mm 17° 161192 - 2,5mm 161193 - 3,5mm 30° 033193 - 3mm 033194 - 4mm	Φ3,5 Torque 30Ncm 020016 Extension 023020 - 0,8mm 023021 - 1,5mm 023022 - 2,5mm 023023 - 3,5mm 023024 - 4,5mm 020129 17° 023122 - 2,5mm 023121 - 2,5mm invert	Φ3,8 Torque 30Ncm 056537 1,5mm 056540 Zero W/ Hex 056542 Zero w/o Hex		Φ3,8 Torque 30Ncm Φ4,5 156104 4.7 3.5 156104 4.7 3.5 306507 - 0,8mm 306508 - 1,5mm 306509 - 2,5mm 306510 - 3,5mm 306511 - 4,5mm 306745 - 0,8x6mm 306500 - 0,8mm 306501 - 1,5mm 306502 - 2,5mm 306503 - 3,5mm 306504 - 4,5mm 306738 - 0,8x6mm			
230009 - Square (Open frame) 230007 Conical (Closed frame) 			103000 Conical (Closed frame) 094000 Square (Open frame) 		100000 - Conical (Closed frame) 023000 - Square (Open frame w/o Hex) 167000 - Square (Open frame W/ Hex) 	021001 Square (Open frame W/ Hex) 	957274 (Open frame) Transfer Scan TS (Ti) 026103 (Open frame) Φ3,6mm 028103 (Closed frame) Φ3,6mm 028105 (Closed frame) Φ4,6mm 					
 957255 Scan Conect MB			 957253 Scan Conect MU		 957259 c/ Hex Scan Conect EC		 957261 Scan Conect UMD		 957232 Scan Conect NP		 957257 - Small Scan Conect TBS 957258 - Large Scan Conect TBL	
 013411			 013410		 013414		 013412		 013401		 013415 Small 013416 Large	
Thread 1.6 Torque 20Ncm 230003			Thread 1.4 Torque 20Ncm 157004 for Metal 158004 Slot for Metal 020157 Exclusive for Zirconia Struc.	Thread 1.4 Torque 20Ncm 011004		Thread 1.4 Torque 20Ncm 020016		Thread 1.6 Torque 30Ncm 156104				
 230006 Multiple titanium			 159004 - Multiple titanium		 112004 - Titanium w/o Hex multiples 170004 - Titanium W/ Hex unitary		 109006 - Titanium Index		 137104 - Titanium Index Φ4,5 137154 - Titanium Index Φ3,5 137124 Linha Zero Titanium w/o Hex 307025 NP Ti Cement E Rotacional			
 055097 CrCo 230004 Ti-short 230011 Ti Cement EX 230010 Calcinable 230005 Ti-long			144001 - Calcinable * Use slotted screw 055096 - CrCo Chrome-Cobalt 144004 TiLaser 036001 TiCement 144008 Ti Cement EX		 105001 - w/o Hex 106001 - W/ Hex 105004 - w/o Hex 055034 - W/ Hex 055095 - w/o Hex 055094 - W/ Hex		Calcinable W/ Hex 107413 Screwed 107417 Cemented 4mm 107422 Cemented 6mm		Pilar 0° Ring Φ4,0 Torque 30Ncm 049211 - 1,5mm 049212 - 2,5mm 049213 - 3,5mm 049214 - 4,5mm 203200 - Esp. 0° Ring 203000 - Esp. Micro 208000 - Esp. Ring 207000 - Micro Ring 203300 - Capsule 0° Ring 203100 - Micro Capsule 201000 - Ring 0° Ring 202000 - Micro 0° Ring		Barra Clips 200100 - Straight Clip 200200 - Clip 35 graus 200300 - Clip 70 graus 204000 - Acrylic Bar 400204 - Titanium Bar Φ2,0 x 25mm	
 230002			 095000		 098000		 099405					

IMPORTANT:

DIGITAL SCANNING AND INFORMATION FOR THE LABORATORY

Its required to correctly identify and inform the laboratory:

- Scan Conexão used (Scan Conect, Transfer Scan or Scan inL)
- about which Implant (NP 24° - FIT 11,5° - GMF 16° - SL - HE RD - HE SD ...), or
- About which Abutment (TiBase S - TiBase L - MU - MB - EC - UMD - UB ...)

SCAN TYPES:

- SC Scan Conect (choose: on implant or on abutment)
- TS Transfer Scan (choose: on implant or on abutment)
- inL Scan inL (only on abutment)


Incorrect identification will result in errors in the positioning of the digital analogue and the prosthesis in general.

www.conexao.com.br/bibliotecas


Unitary - Multiple

Scan inL
Sirona® inLab

Micro Unit
Mini Pilar



Multi Base



TiBase S



TiBase L



DentalCAD

Conexão
Digital Library

